



MEDI CAPTURE

## MediCap<sup>®</sup> MVR Lite

The MediCap<sup>®</sup> MVR Lite is the medical industry's most advanced and compact Full HD Medical Video Recorder. The MVR Lite is the ideal solution to the industry's demand for medical grade recorders that are affordable, robust, and offer the highest quality in video recordings.



HD

1080p  
60



### HD Medical Video Recorder

Record Full HD video from your endoscope, arthroscope, surgical camera, ultrasound, etc. New state-of-the-art recording platform records Full 1080p at 60 frames per second.



### Compact and robust

This rugged MVR features a stainless steel case and fluid resistant front panel. The compact size allows limitless mounting and installation options.



### Intuitive Use

It's so easy, no training is needed. Just plug it in and go. You'll be recording HD videos within minutes of opening the box.

Medical Imaging Made Easy<sup>™</sup>

# Simple Features



## Medical Grade

Designed, manufactured and certified for use in the OR / medical environment.



## Photos

Capture and store HD photos during your video recording.



## Remote control

The included remote allows surgeons to record video, stills, and power on/ off at the touch of a button.



## Foot Pedal

Capture high quality HD video and HD photos to a USB drive using an endoscope trigger or foot pedal (optional).

## Also available...



### MediCap® USB300

All in one full HD/3D recorder with built in monitor, network integration and patient information features

## MediCap® MVR Lite Specs<sup>1</sup>

<b>Device Type</b>	<b>HD Medical Video Recorder</b> Records high-definition video and still images from medical video sources to external USB drive.
<b>Case/Buttons</b>	Stainless steel, fluid resistant
<b>User interface</b>	Overlay on external monitor (record, pause, storage information)
<b>Video Inputs</b>	HDMI, DVI (with adapter)
<b>Video Outputs</b>	HDMI, DVI (with adapter)
<b>Video Formats</b>	H.264 / MPEG4 (1080p60 / 1080p50)
<b>Input video resolutions</b>	1080p50/60, 1080i50/60, 720p50/60, 480p60, 576p50, 1280x1024, 1920x1200p50/60
<b>Image Formats</b>	JPG
<b>External Storage</b>	USB flash drive, USB hard drive
<b>USB Support</b>	USB3.0, high speed, 1 front connector, 2 rear connectors
<b>Footswitch/Trigger</b>	3.5mm mini jack. Compatible with standard footswitches (not included) and closed-contact camera triggers.
<b>Mechanical</b>	Size: 160 x 140 x 40mm / 6.3 x 5.5 x 1.6" Weight: 0.76kg / 1lb 11oz
<b>Environmental</b>	Storage: -40° to +85° C Operating: -20° to +35° C
<b>Power</b>	AC 110-240V, 50/60Hz, 0.4A, 3-prong AC jack
<b>Certifications</b>	CE 93/42/EEC, IEC 60601-1:2005, ANSI/AAMI ES60601-1:2005, CAN/CSA C22.2 No.60601, 1:2008, E361228, FCC Part 15, FDA Class 1, RoHS
<b>Kit Contents</b>	MediCap MVR Lite, USB flash drive, power cable, DVI-to-HDMI adapter, remote control, user guide



<sup>1</sup> Subject to change



Back Panel Connectors